

RNAseq: mRNA to Illumina library

CURRENT STATUS: POSTED

Marcia Sanders
ORFeome

✉ marcia_sanders@unc.edu *Corresponding Author*

Anne Beall
ORFeome

✉ aebell@email.unc.edu *Corresponding Author*

DOI:

10.1038/protex.2015.013

SUBJECT AREAS

Biotechnology *Microbiology*

KEYWORDS

RNAseq, TruSeq Stranded mRNA Sample Preparation, Illumina mRNA library

Abstract

This protocol details the RNAseq sample preparation starting with mRNA and ending with quality control of finished libraries. This protocol was used specifically for virus mRNA samples related to ORFeome lab group.

Introduction

This protocol details the RNAseq sample preparation starting with mRNA. This is a modification of the TruSeq Stranded mRNA Sample Preparation since the published protocol usually starts with Total RNA instead of mRNA. Protocol also includes assays used for quality control of samples.

Reagents

Agilent RNA 6000 Nano Kit (5067-1511)

Agilent High Sensitivity DNA Kit (5067-4626)

Illumina Tru Seq Stranded mRNA Kit Set A (RS-122-2101)

Quant-iT™ PicoGreen® dsDNA Assay Kit (P7589)

Equipment

2100 Electrophoresis Bioanalyzer Instrument (G2939AA)

Fluorescence plate reader

Procedure

1. **Assess mRNA using Agilent RNA 6000 Nano kit**, assay mRNA Nano (Supplementary 1).
 - a. Best if each mRNA sample run in duplicate so an average mRNA concentration can be used.
 - b. Following procedure requires 10-400ng of mRNA for each sample.
2. **Proceed with Illumina TruSeq Stranded mRNA Sample Preparation Guide, Low Sample (LS) Protocol (15031047 Revision E)**.
 - a. Since starting with mRNA, start at page 20 step #12 by adding mix to desired amount of mRNA.
 - i. Illumina recommends having mRNA concentrated into 5 μ L or less to start; however, we have had success with using up to 42 μ L, equal to 400ng of mRNA.
 - b. Proceed with incubation.
 - c. Since starting with mRNA, without beads, skip page 22 step #2. At step #3 only transfer 17 μ L of

mRNA/mix and continue protocol as Illumina has written.

3. **Validate library after complete.** a. Assess quantity of dsDNA library using Quant-iT™ PicoGreen® dsDNA Assay Kit (Supplementary 2). Minimum quantity will depend on laboratory/core the library will be submitted to.
- b. Assess size of library using Agilent High Sensitivity DNA Kit (Supplementary 3) . Refer to the Quick Guide. The final product should be approximately 260bp.

Timing

Agilent RNA 6000 Nano-1.5hrs

TruSeq Stranded mRNA Prep-9hrs

Agilent HS DNA-1.5hrs

Picogreen-0.5hrs

Troubleshooting

Please refer troubleshooting to applicable commercial company technical support.

Anticipated Results

A valid Illumina Library to send for sequencing.

References

1. Agilent High Sensitivity DNA Kit Quick Start Guide-Rev.C. (2013) Available at:
http://www.chem.agilent.com/library/usermanuals/Public/G2938-90322_HighSensitivityDNA_QSG.pdf.
2. Agilent RNA 6000 Nano Kit Quick Start Guide-Rev.C. (2013) Available at:
http://www.chem.agilent.com/library/usermanuals/Public/G2938-90037_RNA6000Nano_QSG.pdf.
3. Low Sample (LS) Protocol. TruSeq Stranded mRNA Sample Preparation Guide. [Revision E] [20-43] (Illumina, Inc. 2013)
4. Quant-iT™ PicoGreen ® dsDNA Reagent and Kits-Rev.10-June-2008. (2008) Available at: <http://tools.lifetechnologies.com/content/sfs/manuals/mp07581.pdf>.

Acknowledgements

UNC High-Throughput Sequencing Facility for sequencing the finished libraries.

Supplementary Files

This is a list of supplementary files associated with this preprint. Click to download.

[Manual_G2938-90037_RNA6000Nano_QSG.pdf](#)

[Manual_G2938-90322_HighSensitivityDNA_QSG.pdf](#)

[Quant-iT_Picogreen_dsDNA.pdf](#)