

MSM Trace Element (TE) solution

Stephan Christel (✉ stephan.christel@lnu.se)

Systems Biology of Extreme Microorganisms, Dopson, Linnaeus University

Mark Dopson


Systems Biology of Extreme Microorganisms, Dopson, Linnaeus University

Method Article

Keywords: mineral, salt, medium, salts, cell, culture, acidophile, bioleaching, biomining, pH, acidophilic, extreme, mineral, medium, media, trace, elements

Posted Date: March 24th, 2015

DOI: <https://doi.org/10.1038/protex.2015.032>

License:   This work is licensed under a Creative Commons Attribution 4.0 International License. [Read Full License](#)

Abstract

Trace element solution for Mineral Salts Medium (MSM, see separate protocol

"<http://www.nature.com/protocolexchange/protocols/3809>":<http://www.nature.com/protocolexchange/protocols/3809>).

To be added when medium is not supplied with mineral. Add 1 ml of 100x Stock TE solution into 99 ml of MSM.

Reagents

Trace Element solution: Ferric chloride hexahydrate, $\text{FeCl}_3 \cdot 6\text{H}_2\text{O}$ Copper sulfate pentahydrate, $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$ Boric acid, HBO_3 Manganese sulfate hydrate, $\text{MnSO}_4 \cdot \text{H}_2\text{O}$ Sodium molybdate dihydrate, $\text{Na}_2\text{MoO}_4 \cdot 2\text{H}_2\text{O}$ Cobalt chloride hexahydrate, $\text{CoCl}_2 \cdot 6\text{H}_2\text{O}$ Zinc sulfate heptahydrate, $\text{ZnSO}_4 \cdot 7\text{H}_2\text{O}$ Sodium selenate, Na_2SeO_4 General: Sulfuric acid Sodium hydroxide

Equipment

Stirring rods, Magnetic stirrer, Measuring colves, pH meter

Procedure

For 100x Stock TE solution: 1100 mg/l $\text{FeCl}_3 \cdot 6\text{H}_2\text{O}$ 50 mg/l $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$ 200 mg/l HBO_3 200 mg/l $\text{MnSO}_4 \cdot \text{H}_2\text{O}$ 80 mg/l $\text{Na}_2\text{MoO}_4 \cdot 2\text{H}_2\text{O}$ 60 mg/l $\text{CoCl}_2 \cdot 6\text{H}_2\text{O}$ 90 mg/l $\text{ZnSO}_4 \cdot 7\text{H}_2\text{O}$ 10 mg/l Na_2SeO_4 - Add the constituents into a final volume of 1 liter of water (pH pre-adjusted to around 2.5), one at a time, dissolving before adding the next. - Adjust the pH to 1.5 with sulfuric acid (and/or sodium hydroxide) -sterile filter

Timing

about 30 minutes

Troubleshooting

-TE solution should be clear and almost colorless. Strong coloration indicates too high pH. -When encountering difficulties during dissolution of a reagent, try adding some more sulfuric acid.