

# Standard Protocols for Flow Cytometry for monocytes, macrophages, DC and T cells

Paloma Riquelme (✉ [paloma.riquelme@klinik.uni-regensburg.de](mailto:paloma.riquelme@klinik.uni-regensburg.de))

Department of Surgery, University Hospital Regensburg

James A. Hutchinson (✉ [james.hutchinson@klinik.uni-regensburg.de](mailto:james.hutchinson@klinik.uni-regensburg.de))

Department of Surgery, University Hospital Regensburg

---

## Method Article

**Keywords:** flow cytometry, Mreg, miTreg, macrophage, T cell, IDO, Foxp3, Helios

**Posted Date:** July 11th, 2018

**DOI:** <https://doi.org/10.1038/protex.2018.064>

**License:**   This work is licensed under a Creative Commons Attribution 4.0 International License.

[Read Full License](#)

---

# Abstract

This document describes standard research protocols used for analysis of marker expression by flow cytometry in the Hutchinson work-group. In Riquelme et al., these protocols have been applied to the analysis of monocyte-derived macrophages, DC and cultured T cells.

# Introduction

This document describes standard research protocols used for analysis of marker expression by flow cytometry in the Hutchinson work-group. In Riquelme et al., these protocols were applied to the analysis of monocyte-derived macrophages and DC, and cultured T cells.

# Reagents

**\*\*General reagents\*\*** Biocoll Separating Solution \ (L6115, Millipore-Biochrom) DPBS without  $\text{Ca}^{2+}$  or  $\text{Mg}^{2+}$  \ (D8537, Sigma) BSA \ (01400.1, Biomol)  $\text{NaN}_3$  \ (S8032, Sigma) EDTA 0.5M \ (15575-038, Gibco Life) Human FcR blocking reagent \ (130-059-901, Miltenyi) 7-AAD \ (559925, BD) Intracellular fixation buffer \ (00-8222, eBiosciences) Permeabilization Buffer \ (00-8333, eBiosciences) FOXP3 Fixation/Permeabilization \ (00-5521-00, eBiosciences) Fixable Live-Dead Viability DyeFluor660 \ (65-0864, eBiosciences) Fixable Live-Dead Viability DyeFluor506 \ (65-0866, eBiosciences) Faser Kit APC \ (130-091-762, Miltenyi) **\*\*Antibodies\*\*** Please see Table 1

# Equipment

FACSCanto II cytometer: Standard configuration with 8 colour, 3 lasers \ (BD Biosciences) DivaTM Software version v.6.1.3 \ (BD Biosciences) Navios cytometer: Standard configuration with 10 colours, 3 lasers \ (Beckman Coulter) Navios clinical software \ (Beckman Coulter) FlowJoTM analysis software v.7.6.5 \ (FlowJo, LLC.) Kaluza analysis software v.1.0 \ (Beckman Coulter)

# Procedure

**\*\*1. REAGENT PREPARATION\*\*** **\*\*1.1 Preparation of FACS buffer\*\*** 1. Make stock solution of 10%  $\text{NaN}_3$  in DPBS 2. Transfer 497ml DPBS to a Schott bottle 3. Add 1 ml 10%  $\text{NaN}_3$  solution to DPBS for an end concentration of 0.02%  $\text{NaN}_3$  4. Add 2 ml 0.5 M EDTA to give an end concentration of 2 mM 5. Add 5 g BSA and allow to dissolve 6. Store solution at 4°C **\*\*2. PROTOCOLS FOR CELL-SURFACE STAINING OF MACROPHAGES, DC AND T CELLS\*\*** **\*\*2.1 Standard extracellular staining protocol\*\*** 1. Harvest cells into cold DBPS.  $2 \times 10^5$  cells are required for each staining reaction. 2. Pellet cells by centrifugation at 300 g for 6 min 3. Aspirate supernatant using vacuum pump 4. Resuspend cells in 100  $\mu\text{l}$ /reaction of FACS buffer + 10% FcR block 5. Incubate at 4°C for 30 min 6. Transfer blocked cell suspensions into FACS tubes at 100  $\mu\text{l}$ /tube 7. Add primary antibodies as indicated 8. Incubate at 4°C in dark for 20 min 9. Add 5  $\mu\text{l}$  7-AAD per reaction. Incubate at 4°C in dark for 10 min 10. Add 2 ml DBPS per tube. Vortex 11. Pellet

cells by centrifugation at 300 g for 6 min 12. Aspirate supernatant using vacuum pump. 13. Resuspend cells in 200 µl DPBS for analysis by flow cytometry. **\*\*2.2 Amplification of PE or APC signals using FASER kit\*\*** 1. Harvest cells into cold DBPS.  $2 \times 10^5$  cells are required for each staining reaction. 2. Pellet cells by centrifugation at 300 g for 6 min 3. Aspirate supernatant using vacuum pump 4. Resuspend cells in 100 µl / reaction of FACS buffer + 10% FcR block 5. Incubate at 4°C for 30 min 6. Transfer blocked cell suspensions into FACS tubes at 100 µl / tube 7. Add primary PE or APC conjugated antibody as indicated 8. Incubate at 4°C in dark for 20 min 9. Add 5 µl 7-AAD per reaction. Incubate at 4°C in dark for 10 min 10. Add 2 ml DBPS per tube. Vortex 11. Pellet cells by centrifugation at 300 g for 6 min 12. Aspirate supernatant using vacuum pump. 13. Resuspend cells in 80 µl FACS buffer + 20 µl FcR block + 10 µl FASER reagent 1 14. Incubate at 4°C in dark for 10 min 15. Wash one time in 2 ml DPBS by centrifugation at 300 g for 6 min 16. Resuspend cells in 80 µl FACS buffer + 20 µl FcR block + 10 µl FASER reagent 2 17. Incubate at 4°C in dark for 10 min 18. Wash one time in 2 ml DPBS by centrifugation at 300 g for 6 min 19. Resuspend cells in 80 µl FACS buffer + 20 µl FcR block + 10 µl FASER reagent 1 20. Incubate at 4°C in dark for 10 min 21. Wash one time in 2 ml DPBS by centrifugation at 300 g for 6 min 22. Resuspend cells in 80 µl FACS buffer + 20 µl FcR block + 10 µl FASER reagent 2 23. Incubate at 4°C in dark for 10 min 24. Wash one time in 2 ml DPBS by centrifugation at 300 g for 6 min 25. Resuspend cells in 200 µl DPBS for analysis by flow cytometry. **\*\*3. INTRACELLULAR STAINING OF HUMAN MREGS FOR IDO\*\*** 1. Harvest cells into cold DBPS.  $5 \times 10^5$  Mregs are required for each staining reaction 2. Pellet cells by centrifugation at 300 g for 6 min 3. Aspirate supernatant using vacuum pump 4. Resuspend cells in 3 ml DBPS 5. Add 0.75 µl Fixable Live-Dead dye 6. Incubate at 4°C in dark for 30 min 7. Wash one time in 10 ml DPBS by centrifugation at 300 g for 6 min 8. Resuspend cells in 300 µl FACS buffer + 10% FcR block for 15 min 9. Transfer blocked cell suspensions into FACS tubes at 100 µl / tube 10. Add primary PE-conjugated CD33 or isotype antibody as indicated in Figure 1 \ (3 reactions) 11. Incubate at 4°C in dark for 30 min 12. Wash samples one time in 2 ml DPBS by centrifugation at 300 g for 6 min 13. Add 100 µl IC Fixation Buffer for 30 min RT in the dark 14. Wash samples twice in 2 ml Perm Buffer by centrifugation at 300 g for 6 min 15. Resuspend samples in 100 µl of FACS buffer + 10% FcR block for 15 min at 4°C 16. Add intracellular staining antibodies as indicated for 30 min at RT in the dark 17. Wash samples once in 2 ml Perm Buffer by centrifugation at 300 g for 6 min 18. Resuspend cells in 200 µl DPBS for analysis by flow cytometry [See figure in Figures section.](#) **\*\*4. INTRACELLULAR STAINING OF HUMAN TREGS FOR FOXP3 AND HELIOS\*\*** **\*\*4.1 Intracellular staining of cultured T cells\*\*** 1. Harvest cells into cold DBPS.  $5 \times 10^5$  cells are required for each staining reaction 2. Pellet cells by centrifugation at 300 g for 6 min 3. Aspirate supernatant using vacuum pump 4. Resuspend cells in 1 ml DBPS 5. Add 1 µl Fixable Live-Dead dye 6. Incubate at 4°C in dark for 30 min 7. Wash one time in 10 ml DPBS by centrifugation at 300 g for 6 min 8. Resuspend cells in 100 µl FACS buffer + 10% FcR block for 15 min 9. Transfer blocked cell suspensions into FACS tubes at 100 µl / tube 10. Add cell surface antibodies as indicated in Figure 2 11. Incubate at 4°C in dark for 30 min 12. Wash samples one time in 2 ml DPBS by centrifugation at 300 g for 6 min 13. Aspirate supernatant using vacuum pump. Gentle vortex 14. Add 500 µl Fixation/Permeabilization 15. Incubate at 4°C in dark for 30 min 16. Wash samples in 2 ml Perm Buffer by centrifugation at 300 g for 6 min 17. Add intracellular staining antibodies \ (FoxP3-PE, optional Helios-APC) as indicated 18. Incubate for 30 min at 4°C in the dark 19. Wash samples 3 times in 2 ml Perm Buffer by centrifugation at 300 g for

6 min 20. Resuspend cells in 200  $\mu$ l Perm Buffer for analysis by flow cytometry [See figure in Figures section](#). **\*\*4.2 Intracellular staining for FoxP3 in T cells from clinical blood samples\*\*** 1 1 ml blood is used for each staining reaction. Dilute blood 1:1 with DBPS 2 Layer blood carefully onto 6 ml Biocoll in a 15 ml tube 3 Centrifuge at 2000 rpm, 20°C, for 25 min with no brake 4 Recover cell layer in the interface 5 Wash with 14 ml DPSB by centrifugation at 1600 rpm for 10 min 6 Aspirate supernatant using vacuum pump 7 Resuspend cells in 100  $\mu$ l / reaction of FACS buffer + 10% FcR block 8 Incubate at 4°C for 30 min 9 Transfer blocked cell suspensions into FACS tubes at 100  $\mu$ l / tube 10 Add cell surface antibodies as indicated in Figure 3 11 Incubate at 4°C in dark for 30 min 12 Add 2 ml DBPS per tube. Vortex 13 Pellet cells by centrifugation at 300 g for 6 min 14 Aspirate supernatant using vacuum pump. Gentle vortex 15 Add 500  $\mu$ l Fixation/Permeabilization 16 Incubate at 4°C in dark for 30 min 17 Wash samples in 2 ml Perm Buffer by centrifugation at 300 g for 6 min 18 Add intracellular staining antibodies as indicated 19 Incubate for 30 min at RT in the dark 20 Wash samples 3 times in 2 ml Perm Buffer by centrifugation at 300 g for 6 min 21 Resuspend cells in 200  $\mu$ l Perm Buffer for analysis by flow cytometry [See figure in Figures section](#). **\*\*5. Gating strategy for identifying TIGIT<sup>+</sup> Tregs among NSG mouse splenocytes\*\*** [See figure in Figures section](#).

## References

1. Hutchinson J.A. et al. MITAP-compliant characterization of human regulatory macrophages. *Transpl Int.* 30(8):765-775 (2017).
2. Hutchinson J.A. et al. Cutting Edge: Immunological consequences and trafficking of human regulatory macrophages administered to renal transplant recipients. *J Immunol.* 187(5):2072-8 (2011)
3. Hutchinson J.A. et al. Human regulatory macrophages. *Methods Mol Biol.* 677:181-92 (2011)

## Figures

		488 Excitation		633 Excitation
		FL2 (PE)	FL4 (PerCP-eFluor710)	FL6 (APC)
PO01	ANTIGEN	CD33	IDO	Live-Dead
	Clone name	P67.6	eyedio	
	Isotype	mIgG1	mIgG1	
	Amount (µl)	5,00	5,00	0,25 ul/ml
	Supplier	BD	eBiosciences	eBio
	Cat. #	347787	46-9477-42	65-0864-14
PO02	ANTIGEN	CD33	Isotype	Live-Dead
	Clone name	P67.6	P3.6.2.8.1	
	Isotype	mIgG1	mIgG1	
	Amount (µl)	5,00	0,30	0,25 ul/ml
	Supplier	BD	eBiosciences	eBio
	Cat. #	347787	9046-4714-120	65-0864-14
PO03	ANTIGEN	Isotype	Isotype	Live-Dead
	Clone name	MOPC-21	P3.6.2.8.1	
	Isotype	mIgG1	mIgG1	
	Amount (µl)	10,00	0,30	0,25 ul/ml
	Supplier	BD	eBiosciences	eBio
	Cat. #	550854	9046-4714-120	65-0864-14

Figure 1

Staining panel IDO

		488 Excitation				633 Excitation		405 Excitation		
		FL1 (FITC)	FL2 (PE)	FL3 (PC5.5)	FL4 (PE-Cy7)	FL5 (APC)	FL6 (APC-H7)	FL7 (V450)	FL8 (Aqua)	
TREG_PANEL	ANTIGEN	CD8	FOXP3	CD3	TIGIT	Drop-in	CD25	CD4	Live-Dead	
	Clone name	SK1	PCH101	SK7	MBSA43		M-A251	RPA-T4		
	Isotype	mIgG1k	rat IgG2ak	mIgG1k	mIgG1		mIgG1k	mIgG1k		
	Amount (µl)	20 ul	10 ul	20 ul	10 ul		5 ul	5 ul	+	
	Supplier	BD	eBio	BD	eBio		BD	BD	eBio	
		Cat. #	345772	12-4776-42	332771	25-9500-42		560225	560345	
		Status	CE/IVD	RUO	CE/IVD	RUO		RUO	RUO	


Figure 2

Staining panel iTreg

		488 Excitation				633 Excitation		405 Excitation	
		FL1 (FITC)	FL2 (PE)	FL3 (PC5.5)	FL4 (PE-Cy7)	FL5 (APC)	FL6 (APC-H7)	FL7 (V450)	FL8 (Aqua)
TREG_01	ANTIGEN	CD8	FOXP3	CD3	TIGIT	CD226	CD25	CD4	CD45
	Clone name	SK1	PCH101	SK7	MBSA43	DNAX11	M-A251	RPA-T4	2D1
	Isotype	mlgG1k	rat IgG2ak	mlgG1k	mlgG1	mlgG1k	mlgG1k	mlgG1k	mlgG1k
	Amount (µl)	20 ul	10 ul	20 ul	10 ul	10 ul	5 ul	5 ul	5 ul
	Supplier	BD	eBio	BD	eBio	Miltenyi	BD	BD	BD
	Cat. #	345772	12-4776-42	332771	25-9500-42	130-092-477	560225	560345	655873
Status	CE/IVD	RUO	CE/IVD	RUO	RUO	RUO	RUO	CE/IVD	
TREG_02	ANTIGEN	CD8	FOXP3	CD3	TIGIT	Isotype	CD25	CD4	CD45
	Clone name	SK1	PCH101	SK7	MBSA43	IS5-21F5	M-A251	RPA-T4	2D1
	Isotype	mlgG1k	rat IgG2ak	mlgG1k	mlgG1	mlgG1k	mlgG1k	mlgG1k	mlgG1k
	Amount (µl)	20 ul	10 ul	20 ul	10 ul	10 ul	5 ul	5 ul	5 ul
	Supplier	BD	eBio	BD	eBio	Miltenyi	BD	BD	BD
	Cat. #	345772	12-4776-42	332771	25-9500-42	130-092-214	560225	560345	655873
Status	CE/IVD	RUO	CE/IVD	RUO	RUO	RUO	RUO	CE/IVD	
TREG_03	ANTIGEN	CD8	FOXP3	CD3	TIGIT	CD278	CD25	CD4	CD45
	Clone name	SK1	PCH101	SK7	MBSA43	ISA-3	M-A251	RPA-T4	2D1
	Isotype	mlgG1k	rat IgG2ak	mlgG1k	mlgG1	mlgG1k	mlgG1k	mlgG1k	mlgG1k
	Amount (µl)	20 ul	10 ul	20 ul	10 ul	5 ul	5 ul	5 ul	5 ul
	Supplier	BD	eBio	BD	eBio	eBio	BD	BD	BD
	Cat. #	345772	12-4776-42	332771	25-9500-42	17-9948-42	560225	560345	655873
Status	CE/IVD	RUO	CE/IVD	RUO	RUO	RUO	RUO	CE/IVD	

Figure 3

Staining panel iTreg in peripheral blood



**Figure 4**

Gating strategy for identifying TIGIT<sup>+</sup> Tregs among NSG mouse splenocytes

## Supplementary Files

This is a list of supplementary files associated with this preprint. Click to download.

- [supplement0.pdf](#)