

# Standard protocols for immune profiling of peripheral blood leucocyte subsets by flow cytometry using DuraClone IM reagents

**Katharina Kronenberg**

Department of Surgery, University Hospital Regensburg, Germany

**Paloma Riquelme** (✉ [Paloma.Riquelme@ukr.de](mailto:Paloma.Riquelme@ukr.de))

Department of Surgery, University Hospital Regensburg, Germany <https://orcid.org/0000-0001-8962-3156>

**James A. Hutchinson**

Department of Surgery, University Hospital Regensburg, Germany <https://orcid.org/0000-0002-1199-4210>

---

## Method Article

**Keywords:** Flow Cytometry, Immune monitoring, DuraClone, IM Phenotyping Basic Tube, IM T cell subsets Tube, IM TCRs Tube, IM Granulocytes Tube, IM Dendritic Cells Tube, IM B cells Tube, IM Treg Tube

**Posted Date:** March 4th, 2021

**DOI:** <https://doi.org/10.21203/rs.3.pex-757/v1>

**License:**  This work is licensed under a Creative Commons Attribution 4.0 International License.

[Read Full License](#)

---

# Abstract

This document describes standard protocols used by the Hutchinson group for profiling human peripheral blood leucocytes by flow cytometry analyses using DuraClone IM reagents.

## Introduction

This document describes standard protocols used by the Hutchinson group for analysis of human blood leucocyte populations by flow cytometry (1). DuraClone reagents are dry pre-formulated antibody panels for extensive phenotyping of T cells, B cells, NK cells, dendritic cells and granulocytes in whole blood. In our experience, DuraClone IM Tubes are a convenient, reliable and highly standardised option for immune monitoring of clinical trials. In part, these protocols were adapted from the publications of Streitz, Kverenland and colleagues (2-3), as well as manufacturer's recommendations issued by Beckman Coulter.

## Reagents

### Reagents

DuraClone IM Phenotyping Basic Tube (B53309, Beckman Coulter)

DuraClone IM T cell subsets Tube (B53328, Beckman Coulter)

DuraClone IM TCRs Tube (B53340, Beckman Coulter)

DuraClone IM Granulocytes Tube (B88651, Beckman Coulter)

DuraClone IM Dendritic Cells Tube (B53351, Beckman Coulter)

DuraClone IM B cells Tube (B53318, Beckman Coulter)

DuraClone IM Treg Tube (B53346, Beckman Coulter)

VersaLyse Solution (A09777, Beckman Coulter)

DPBS without  $\text{Ca}^{2+}$  or  $\text{Mg}^{2+}$  (D8537, Sigma)

IOtest 3 Fixative Solution (10x) (A07800, Beckman Coulter)

FBS South American (10270106, Thermo Fisher), heat-inactivated

Ecotainer Aqua B, Distilled water (0082423E, Braun)

PerFix-nc Kit (B31168, Beckman Coulter)

## Materials

S-Monovette® 2,7 ml, K3 EDTA (05.1167, Sarstedt)

15 ml tube, PP, 17/120 mm, conical bottom (188271, Greiner Bio-One)

## Equipment

### Equipment

Navios cytometer: Standard configuration with 10 colours, 3 lasers (Beckman Coulter)

Cytometry List Mode Data Acquisition and Analysis Software (Beckman Coulter)

Kaluza analysis software v.1.3 (Beckman Coulter)

## Procedure

### 1. Preparation of reagents

#### 1.1 Preparation of 0.8% IOTest 3 Fixative Solution

1. Dilute 100 µl IOTest 3 Fixative Solution (10X) with 900 µl PBS to obtain 0.8 % IOTest 3 Fixative Solution

#### 1.2 Preparation of 0.1% IOTest 3 Fixative Solution

1. Dilute 12.5 µl IOTest 3 Fixative Solution (10X) with 987.5 µl PBS to obtain 0.1 % IOTest 3 Fixative Solution

#### 1.3 Preparation of 1X PerFix-nc Buffer 3

1. Prepare a 1/10 dilution of 10X Buffer 3 in ddH<sub>2</sub>O: dilute 1 part 10X Buffer 3 with 9 parts of ddH<sub>2</sub>O

### 2. Sample preparation using DuraClone IM Tubes: Phenotyping Basic, T cell subsets, TCRs and Dendritic Cells

1. Collect peripheral blood samples directly into EDTA-vacutainers
2. Samples can be stored at 4 °C for a maximum of 4 h before processing
3. Mix the blood sample by repeatedly inverting the vacutainer
- 4.
- 4.1. For DuraClone IM Phenotyping Basic, T cell, TCRs and Granulocytes Tube: transfer 100 µl blood into the DuraClone Tube
- 4.2. For DuraClone IM Dendritic Cells Tube: transfer 200 µl blood into the DuraClone Tube
5. Vortex tube for 8 seconds
6. Incubate at room temperature in the dark for 15 min
7. Add 2 ml VersaLyse Solution
8. Vortex tube for 3 seconds
9. Incubate at room temperature in the dark for 15 min
10. Pellet cells by centrifugation at 200 g for 5 min
11. Aspirate supernatant using a vacuum pump
12. Vortex tube for 3 seconds
13. Add 3 ml DPBS
14. Pellet cells by centrifugation at 200 g for 5 min
15. Aspirate supernatant using a vacuum pump
16. Resuspend cell pellet in 380 µl 0.8 % IOTest 3 Fixative Solution for data collection using a Navios flow cytometer

### **3. Sample preparation using DuraClone IM Tubes: B cells**

1. Collect peripheral blood samples directly into EDTA-vacutainers
2. Samples can be stored at 4 °C for a maximum of 4 h before processing
3. Mix the blood sample by repeatedly inverting the vacutainer

4. Transfer 10 ml DPBS into a 15 ml conical bottom tube
5. Add 300  $\mu$ l blood to the 10 ml DPBS
6. Mix by inverting the tube
7. Pellet cells by centrifugation at 300 g for 10 min
8. Aspirate supernatant using a vacuum pump
9. Vortex tube for 3 seconds
10. Add 10 ml DPBS
11. Mix by inverting the tube
12. Pellet cells by centrifugation at 300 g for 5 min
13. Aspirate supernatant using a vacuum pump
14. Resuspend cell pellet in DPBS to get a final volume of 300  $\mu$ l
15. Transfer 100  $\mu$ l of washed blood sample into the DuraClone IM B cells Tube
16. Vortex tube for 8 seconds
17. Incubate at room temperature in the dark for 15 min
18. Add 2 ml VersaLyse
19. Vortex for 3 seconds
20. Incubate at room temperature in the dark for 15 min
21. Pellet cells by centrifugation at 200 g for 5 min
22. Aspirate supernatant using a vacuum pump
23. Vortex tube for 3 seconds
24. Add 3 ml DPBS
25. Pellet cells by centrifugation at 200 g for 5 min
26. Aspirate supernatant using vacuum a pump
27. Resuspend cells in 380  $\mu$ l 0.1 % IOTest 3 Fixative Solution for data collection using a Navios flow cytometer

#### 4. Sample preparation using DuraClone IM Tubes: Tregs

1. Collect peripheral blood samples directly into EDTA-vacutainers
2. Samples can be stored at 4 °C for a maximum of 4 h before processing
3. Mix the blood sample by repeatedly inverting the vacutainer
4. Transfer 50 µl blood into the DuraClone IM Treg Tube 1
5. Vortex for 8 seconds
6. Incubate at room temperature in the dark for 15 min
7. Add 3 ml DPBS
8. Pellet cells by centrifugation at 500 g for 5 min
9. Aspirate supernatant using vacuum pump
10. Vortex for 3 seconds
11. Add 50 µl heat-inactivated FCS
12. Mix by pipetting
13. Add 5 µl PerFix-nc Buffer 1
14. Vortex for 8 seconds
15. Incubate at room temperature in the dark for 15 min
16. Add 400 µl PerFix-nc Buffer 2
17. Vortex for 8 seconds
18. Transfer content of DuraClone IM Treg Tube 1 into DuraClone IM Treg Tube 2 by pipetting
19. Vortex Tube 2 for 8 seconds
20. Incubate at room temperature in the dark for 60 min
21. Add 3 ml DPBS
22. Incubate at room temperature in the dark for 5 min

23. Pellet cells by centrifugation at 500 g for 5 min
24. Aspirate supernatant using a vacuum pump
25. Vortex for 3 seconds
26. Add 3 ml 1X PerFix-nc Buffer 3
27. Pellet cells by centrifugation at 500 g for 5 min
28. Aspirate supernatant using a vacuum pump
29. Resuspend cells in 380 µl 1X PerFix-nc Buffer 3 for data collection using a Navios flow cytometer

### **Analysis of data obtained using DuraClone IM Tubes**

Template manual gating strategies for extracting leucocyte subset frequencies from flow cytometry data using Kaluza software are given in Figure 1 to 7. Analyses are ideally performed by experienced operators blinded to clinical outcomes.

Fig. 1 DuraClone IM Phenotyping Basic Tube

Fig. 2 DuraClone IM T cell subsets Tube

Fig. 3 DuraClone IM TCRs Tube

Fig. 4 DuraClone IM Granulocytes Cells Tube

Fig. 5 DuraClone IM Dendritic Cells Tube

Fig. 6 DuraClone IM B cells Tube

Fig. 7 DuraClone IM Treg Tube

## **Troubleshooting**

## **Time Taken**


## **Anticipated Results**

## **References**

1. Hutchinson, JA. Predicting Early Viral Control under Direct-Acting Antiviral Therapy for Chronic Hepatitis C Virus Using Pretreatment Immunological Markers. *Front Immunol.* 7;9:146 (2018)
2. Kverneland, AH. Age and gender leucocytes variances and references values generated using the standardized ONE-Study protocol. *Cytometry A.* 89(6):543-64 (2016)
3. Streitz, M. et al. Standardization of whole blood immune phenotype monitoring for clinical trials: panels and methods from the ONE study. *Transplant Res.* 25;2(1):17 (2013)

## Acknowledgements

## Figures




**Figure 1**

DuraClone IM Phenotyping Basic Tube




**Gating Strategy for the DuraClone IM T cell subsets Tube**



**Figure 2**

DuraClone IM T cell subsets Tube


**Gating Strategy for the DuraClone IM TCRs Tube**



**Figure 3**

DuraClone IM TCRs Tube


**Gating Strategy for the DuraClone IM Granulocytes Tube**



**Figure 4**

DuraClone IM Granulocytes Cells Tube


**Gating Strategy for the DuraClone IM Dendritic Cells Tube**



**Figure 5**

DuraClone IM Dendritic Cells Tube


**Gating Strategy for the DuraClone IM B cells Tube**



**Figure 6**

DuraClone IM B cells Tube

**Gating Strategy for the DuraClone IM Treg Tube**



**Figure 7**

DuraClone IM Treg Tube