

Immunostaining of DDR proteins in senescent cells

Raffaella Di Micco

FIRC Institute of Molecular Oncology

Marzia Fumagalli

FIRC Institute of Molecular Oncology

Fabrizio d'Adda di Fagagna

FIRC Institute of Molecular Oncology

Method Article

Keywords: H2AX, ATMpS1981, pS/TQ, SMC1pS996, immunofluorescence, DNA damage checkpoints, DNA damage response, senescence, oncogene-induced senescence

Posted Date: December 6th, 2006

DOI: <https://doi.org/10.1038/nprot.2006.367>

License:  This work is licensed under a Creative Commons Attribution 4.0 International License.

[Read Full License](#)

Abstract

Introduction

We present here a brief yet detailed methodology to study the activation of a number of DNA damage response factors in senescent human fibroblasts.

Reagents

poly-D-lysine \(\text{SIGMA P6407 5 mg}\) Methanol/Aceton 1:1 PBG 10x stock solution: 5% BSA, 2% of gelatin from cold water fish skin \(\text{SIGMA cat. G-7765}\) in PBS 1x Antibodies: γ H2AX mouse Upstate Biotechnology \(\text{cat.05-636}\) 1:200 ATMpS1981 rabbit Rockland \(\text{cat. 600-401-400}\) 1:300 pS/TQ rabbit Cell Signaling Technology \(\text{cat. \#2851}\) 1:100 SMC1pS996 rabbit Bethyl Laboratories \(\text{cat. A300-050A}\) 1:500 DAPI \(\text{SIGMA cat.D-9564}\) Mowiol mounting medium \(\text{CALBIOCHEM cat. 475904}\)

Procedure

- 1) Prepare poly-D-lysinated coverslips \(\text{poly D lysine } 50\mu\text{g/ml}\)
- 2) Plate $15\text{-}20 \times 10^3$ cells on coverslips 1 day before the staining
- 3) Wash cells gently $2 \times 5'$ in PBS1x
- 4) Fix and permeabilize cells for 2' with methanol/acetone 1:1 at RT
- 5) Block cells for 30' with PBG1x
- 6) Stain with primary antibodies diluted in PBG1x at indicated concentration for 1h at RT in a humidified chamber
- 7) Wash cells gently $3 \times 5'$ in PBG1x
- 8) Stain with secondary antibodies diluted in PBG1x for 1h at RT in a dark humidified chamber
- 9) Wash cells gently $2 \times 5'$ in PBG1x and $2 \times 5'$ in PBS1x
- 10) Stain with Dapi for 5' in the dark
- 11) Mount coverslips with $3\mu\text{l}$ of Mowiol mounting medium
- 12) Let the coverslips dry before microscope analysis

Timing

4 h

Critical Steps

Use always a positive \(\text{X-Rays treated cells } 2\text{-}5 \text{ Gy}\) and a negative control \(\text{proliferating or even better contact inhibited \(\text{(for a short period) quiescent cells}\) to reduce background}\)

Troubleshooting

Using different type of fixative can allow better staining depending on antibody Lipofuscin accumulation in senescent cells can give an high cytosolic autofluorescent background material

References

Figures




Figure 1

Confocal analysis of DDR proteins (γ H2AX=green and 53BP1=red) in oncogene-induced senescent cells. In blue dapi staining.