

Protocol for lentiviral knock down in mouse ES cells

W.W.M. Pim Pijnappel (✉ w.pijnappel@erasmusmc.nl)

Pijnappel/Schöler/Timmers

Marijke P.A. Baltissen

H.T.Marc Timmers

Method Article

Keywords: mouse embryonic stem cells, lentiviral knock down, TFIIID

Posted Date: March 26th, 2013

DOI: <https://doi.org/10.1038/protex.2013.036>

License: © ⓘ This work is licensed under a Creative Commons Attribution 4.0 International License.

[Read Full License](#)

Abstract

This protocol describes the procedure for lentiviral-mediated knock down of transcription factors like TFIIID in mouse embryonic stem cells, and how this can be used to test for factors that are required to maintain the pluripotent state.

Procedure

Note It is important to carefully check the integrity of the shRNA-pLKO DNA for each DNA preparation by restriction enzyme digestion. The reason is that the short hairpin sequence is easily removed from the plasmid by recombination during propagation in DH5alpha cells. Recombination may be avoided by using fresh antibiotics and a limited number of bacterial cell divisions. Day 1 - morning: plate Cos-7 cells in Cos-7 medium at ~40% confluency on 10 cm dishes. 1 dish per transfection. Make sure cells are distributed evenly over the plate. - Afternoon: transfect plated Cos-7 cells as follows: • Tube 1: Mix 300 ul DMEM low glucose (Lonza) with 24 ul Fugene 6 by flicking. Incubate 5 min at RT • Tube 2: Mix the following plasmids: o 2 µg pRSV-rev o 2 µg pMDLg/pRRE o 2 µg pMD2.G o 2 µg pLKO-short hairpin construct • add contents of Tube 1 to Tube 2, mix by flicking, incubate 20-30 min at RT. • Add mixture to Cos-7 cells. Day 2 - replace the medium of the transfected Cos-7 cells with 10 ml ES medium/10 cm dish. - Plate IB10 ES cells on gelatinized dishes at 4.4×10^4 cells/10 cm dish Day 3 - collect viral supernatant from Cos-7 cells using a syringe and filter through 0.45 µm low protein binding filter (to remove cells). Add fresh ES medium to the Cos-7 cells. - Remove the medium from the IB10 cells - Supplement the viral supernatant with 8 µg /ml final concentration polybrene, mix, and add to the IB10 cells. 10 ml undiluted viral supernatant is used to infect 1 x 10 cm dish of IB10 cells. Day 4 - repeat the infection as described for day 3. Add fresh ES medium to the Cos-7 cells. Day 5 - Optional: repeat the infection once more as described for day 3. Discard the Cos-7 cells. 1-2 days after the last infection: - From this point, refresh the medium of the infected IB10 cells every day with ES medium containing 1.8 µg /ml puromycin. Note: it is important to assess the condition of the IB10 cells. They should have clearly grown compared to the day of plating (day 2), and have formed small colonies. If this is not the case 24 hr after the last infection, it is better to wait an additional day before starting the puromycin selection. Day 11 - analyze the cells for alkaline phosphatase colony staining or RT-qPCR to determine knock down efficiency. Medium composition: Cos-7 medium - 500 ml DMEM low glucose (Lonza BE12-707F) - 55 ml FBS (Hyclone cat SV30160.03) - 5 ml Penicillin/Streptomycin (Lonza cat 17-603E) - 5 ml glutamine (Lonza cat 17-605E) ES cell medium - 500 ml DMEM high glucose (Lonza BE 12-604F) - 90 ml FBS (Hyclone cat SV30160.03) - 6 ml Penicillin/Streptomycin (Lonza cat 17-603E) - 6 ml Glutamine (Lonza cat 17-605E) - 6 ml sodium pyruvate (Lonza BE 13-115E) - 6 ml beta-mercaptoethanol (stock: 70 ul of beta mercaptoethanol (Sigma M-6250) in 100 ml H₂O, aliquot and store -20C) - 6 ml non essential amino acids (Lonza 13-114E) - 60 ul LIF (500 U/ml final) (Gibco 13275-029)

Timing

11 days

Anticipated Results

Knock down efficiencies range between 50-80% and should be checked by RT-qPCR analysis. Puromycin selection will ensure a continued knock down of the target mRNA. Knock down of many TFIIID subunits in mouse ES cells will result in colonies with a differentiated morphology and with reduced alkaline phosphatase staining.

References

1. Pijnappel WW, Esch D, Baltissen MP, Wu G, Mischerikow N, Bergsma AJ, van der Wal E, Han DW, Bruch HV, Moritz S, Lijnzaad P, Altelaar AF, Sameith K, Zaehres H, Heck AJ, Holstege FC, Schöler HR, Timmers HT. A central role for TFIIID in the pluripotent transcription circuitry. *Nature*. 2013 Mar 17. doi: 10.1038/nature11970. [Epub ahead of print] PMID: 23503660 [PubMed - as supplied by publisher]

Acknowledgements

We thank M. de Bruijn and O. Kranenburg for advice on lentiviral knockdown experiments