

Isolation of exosomal RNA from serum or plasma using ultracentrifugation and the Qiagen miRNeasy Micro kit

Roopali Gandhi (✉ rgandhi@rics.bwh.harvard.edu)

Department of Neurology, Center for Neurologic Diseases, Brigham and Women's Hospital, Harvard Medical School, Boston, MA, USA

Louise C. Laurent (✉ [llaurent@ucsd.edu](mailto:lolaurent@ucsd.edu))

Department of Reproductive Medicine, University of California, San Diego, La Jolla, CA, USA

Xandra O. Breakefield (✉ breakefield@hms.harvard.edu)

Departments of Neurology and Radiology and Program in Neuroscience, Harvard Medical School, Massachusetts General Hospital, Boston, MA, USA

Leonora Balaj

Department of Neurology and Program in Neuroscience, Harvard Medical School, Massachusetts General Hospital, Boston, MA, USA

Anu Paul

Department of Neurology, Center for Neurologic Diseases, Brigham and Women's Hospital, Harvard Medical School, Boston, MA, USA

Parham Nejad

Department of Neurology, Center for Neurologic Diseases, Brigham and Women's Hospital, Harvard Medical School, Boston, MA, USA

Srimeenakshi Srinivasan

Department of Reproductive Medicine, University of California, San Diego, La Jolla, CA, USA

Xuan Zhang

Department of Neurology, Harvard Medical School, Massachusetts General Hospital, Boston, MA, USA

Method Article

Keywords: exRNA, extracellular vesicles, exosomes, serum, plasma, blood, ultracentrifugation, RNA isolation, miRNeasy kit

Posted Date: June 16th, 2017

DOI: <https://doi.org/10.1038/protex.2017.079>

License:  This work is licensed under a Creative Commons Attribution 4.0 International License.

[Read Full License](#)

Abstract

This protocol describes how to isolate exosomal RNA from serum/plasma by using ultracentrifugation for exosomal enrichment followed by the Qiagen miRNeasy Micro kit for RNA isolation. It supercedes an earlier protocol (<http://dx.doi.org/10.1038/protex.2015.111>). The new protocol is based on our experience vetting the original protocol in multiple labs. The major change is a decrease in input volume of biofluid from 1 mL to 500 μ L in order to match the input volume of other protocols the ERC consortium is developing in tandem.

Introduction

Extracellular RNAs (exRNAs) have been found to play an important role in intercellular communication in the body. exRNAs are present in biofluids, in complexes with lipoproteins and ribonucleoproteins, and in extracellular vesicles such as exosomes and microvesicles. Further study and characterization of exRNAs and their carriers could lead to identification of new biomarkers and have potential for development of novel therapeutics.

Reagents

Qiagen miRNeasy micro kit (217084) Chloroform (Sigma-Aldrich – 319988) 100% ethanol (Koptec – V1016) 70% ethanol RNase-free water (Ambion – AM9937) Phosphate-buffered saline (PBS)

Equipment

Ultracentrifuge Swinging bucket rotor (TLA110/MLA5) Microfuge 1.5 mL Microfuge tubes Vortexer Phase lock gel tubes, 2 mL (VWR – 10847-802)

Procedure

1. Start with 500 μ L serum or plasma.
2. Bring up volume to fill ultracentrifuge tube (2.3-2.5 mL) with PBS.
3. Centrifuge for 70 min at 100,000 x g at 4°C.
4. Discard supernatant and resuspend pellet in PBS to fill ultracentrifuge tube.
5. Centrifuge for 70 min at 100,000 x g at 4°C.
6. Discard supernatant.
7. Resuspend pellet in 50 μ L PBS and store at -80°C, or
8. Proceed to miRNeasy RNA isolation by adding 700 μ L Qiazol to pellet.

Acknowledgements

PN, AP and RG were supported by NIH grant UH2-UH3 TR000890. SS and LL were supported by U01 HL126494 and UH3 TR000906. XZ, LB and XOB were supported by U19 CA179463.