

Protein A/G-based enzyme-linked immunosorbent assay for detection of anti-*Pythium insidiosum* antibodies in human and animal subjects

Chalisa Jaturapaktrarak

Mahidol University Faculty of Medicine Ramathibodi Hospital

Penpan Payattikul

Mahidol University Faculty of Medicine Ramathibodi Hospital

Tassanee Lohnoo

Mahidol University Faculty of Medicine Ramathibodi Hospital

Yothin Kumsang

Mahidol University Faculty of Medicine Ramathibodi Hospital

Aree Laikul

Kasetsart University

Watcharapol Pathomsakulwong

Kasetsart University

Chompoonek Yurayart

Kasetsart University

Walaiporn Tonpitak

Mahanakorn University of Technology

Theerapong Krajaejun (✉ mr_en@hotmail.com)

Mahidol University Faculty of Medicine Ramathibodi Hospital <https://orcid.org/0000-0003-0545-3765>

Research note

Keywords: *Pythium insidiosum*; Pythiosis; Immunodiagnosis; ELISA; Protein A/G

Posted Date: February 17th, 2020

DOI: <https://doi.org/10.21203/rs.2.22143/v2>

License:   This work is licensed under a Creative Commons Attribution 4.0 International License.

[Read Full License](#)

Version of Record: A version of this preprint was published at BMC Research Notes on March 6th, 2020.
See the published version at <https://doi.org/10.1186/s13104-020-04981-y>.

Abstract

Objectives: Pythiosis is a deadly infectious disease caused by *Pythium insidiosum*. Reports of both human and animal pythiosis are on the rise worldwide. Prognosis of the pythiosis patients relies on early diagnosis and prompt treatment. There are needs for an immunodiagnostic test that can detect the disease in both humans and animals. This study aims at reporting an optimized protocol for the development of a protein A/G-based enzyme-linked immunosorbent assay (ELISA) for the detection of anti- *P. insidiosum* antibody in multiple host species. **Results:** A total of 25 pythiosis and 50 control sera, obtained from humans, horses, dogs, cats, and cows, were recruited for the assay development. With a proper ELISA cutoff point, all pythiosis sera can ultimately be distinguished from the control sera. The successfully-developed protein A/G-based ELISA can detect the anti- *P. insidiosum* antibodies in serum samples of both humans and animals. It is a versatile, feasible-to-develop, and functional immunodiagnostic assay for pythiosis.

Introduction

Pythiosis is a deadly infectious disease caused by the pathogenic oomycete *Pythium insidiosum* [1, 2]. The disease affects humans and animals, especially horses and dogs. Geographic distribution of pythiosis covers the tropical and subtropical countries [2]. Pythiosis demonstrates high rates of morbidity and mortality [3]. Most of the affected individuals lose an infected organ, and many patients die from the disease. Prognosis of the pythiosis patients relies on early diagnosis and prompt treatment. Establishing the definitive diagnosis requires a reliable diagnostic method, i.e., organism isolation and identification, immunodiagnostic test, molecular assay, and proteomic analysis [1, 2, 4–13].

Immunodiagnostic tests for pythiosis have gained popularity due to simplicity and short turnaround time. Various immunodiagnostic methods, such as immunodiffusion (ID), enzyme-linked immunosorbent assay (ELISA), hemagglutination (HA), immunochromatographic test (ICT), have been developed to facilitate the diagnosis of pythiosis [10–20]. Each immunodiagnostic test has some advantages/disadvantages over the others. For example, ID and HA detect the antibodies in both humans and animals but possess limited sensitivity. ELISA is a multi-step assay but exhibits relatively-high diagnostic performance. ICT is a rapid and user-friendly assay, but the development of this test is complicated. Because reports of both human and animal pythiosis are on the rise worldwide, there are needs for an immunodiagnostic test that can detect the disease in a broad range of host species. Due to the lack of a versatile, feasible-to-develop, and functional immunodiagnostic assay, the current study aims at reporting an optimized protocol to develop an ELISA for the detection of anti-*P. insidiosum* antibody in humans and other animals. The assay relied on the use of protein A/G, which is a molecule that can bind the immunoglobulins from various animal species [16].

Methods

We extracted proteins from the *P. insidiosum* strain Pi-S [21, 22], using the established method [17, 23]. In short, the organism was maintained on Sabouraud dextrose agar [10 mg/ml peptone (Becton, Dickinson and Company, Maryland, USA), 20 mg/ml dextose (HiMedia Laboratories, Mumbai, India) and 15 mg/ml agar (Becton, Dickinson and Company, Maryland, USA)]. Ten excised agar cubes (~0.5 cm³) with *P. insidiosum* colony were incubated with shaking (150 rpm) at 37 °C for 10 days in 100 ml of Sabouraud dextrose broth. The growing organism was removed from the cultured broth by consecutive filtrations through a filter paper (Whatman No.1, GE healthcare, Belfast, UK) and a 0.22-µm-pore-size Durapore membrane (Merck Millipore, Cork, Ireland). The filtrated broth was concentrated by centrifugation (10,000x g) with an Amicon Ultra 15M tube (Merck Millipore, Cork, Ireland). The concentrated broth, now called culture filtrate antigen (CFA), was stored at -30 °C.

We developed a multiple host-specific ELISA by modifying the human-specific ELISA protocols of Chareonsirisuthigul et al [15] and Lohnoo et al [23]. In brief, a 96-well flat-bottom polystyrene plate (Corning, New York, USA) was coated overnight at 4 °C with 100 µl/well of CFA (5 µg/ml) in 0.1 M carbonate buffer pH 9.6 [0.2 M Na₂CO₃ and 0.2 M NaHCO₃ (Merck, Darmstadt, Germany)] and 1.5% NaCl (Merck, Darmstadt, Germany). The plate was washed 4 times with the washing buffer [phosphate-buffered saline pH 7.4 (PBS); 137 mM NaCl, 2.7 mM KCl, 10 mM Na₂HPO₄, and 1.76 mM KH₂PO (Merck, Darmstadt, Germany) with 0.05% Tween-20 (Calbiochem, California, USA)], blocked with 250 µl of 0.5% bovine serum albumin (Merck, Darmstadt, Germany) in PBS at 37 °C for 60 min, and washed 4 times again. Each serum (1 ul) was added to 9 ul of PBS. Then, the diluted serum (5 ul) was mixed with 795 ul of PBS. The resulting serum sample (1:1,600 in PBS; 100 µl) was added to each well and incubated at 37 °C for 60 min. After washing (as above), 100 µl of protein A/G conjugated with peroxidase (Bio-Rad Laboratories, California, USA) (1:100,000 in PBS) was added to each well and stored at 37 °C for 60 min. The ELISA plate was washed once again. Chromogenic substrate (3,3',5,5'-Tetramethylbenzidine and H₂O₂; Thermo Scientific, Rockford, USA) was applied to each well (100 µl/well) and incubated at room temperature for 3 min in a dark chamber. Sulfuric acid (0.5 N; 100 µl/well) was used to stop the enzymatic reaction. Optical density (OD) of each serum sample was measured at the 450-nm wavelength.

A total of 25 serum samples were derived from 20 Thai patients, 4 horses, and a dog with the clinical features of pythiosis [1–3], together with either successful detection of *P. insidiosum* by culture and molecular analysis [6, 24–26] or serum anti-*P. insidiosum* antibodies by ICT [16]. Fifty control (negative) sera were obtained from 10 each of healthy individuals, horses, dogs, cats, and cows with no clinical sign of pythiosis in the absence of serum anti-*P. insidiosum* antibodies by ICT. Positive (from a known pythiosis patient) and negative (from a healthy blood donor) control sera were analyzed in parallel. All serum samples were tested in duplicate by the protein A/G-based ELISA. ELISA value (EV) was an OD value of a serum sample divided by an OD value of the negative control serum. ELISA cutoff value was derived from the mean EV of the control sera plus 2, 3, 4, or 5 standard deviations (SD). Diagnostic performance (i.e., specificity, sensitivity, and accuracy) were calculated as described previously [15] (**Table 1**). Distribution of EVs (representing a detectable level of anti-*P. insidiosum* antibodies) of the pythiosis (n

= 25) and control (n = 50) sera from both humans and animals was shown in **Figure 1A**. Mean EV of the pythiosis sera was 41.0 (range, 23.1-46.9), whereas that of the control sera was 3.0 (range, 0.4-12.2). The “mean EV+5SDs” cutoff point (EV, 13.4) provided the best sensitivity, specificity, and accuracy for the detection of the anti-*P. insidiosum* antibodies (**Table 1**). Furthermore, protein A/G-based ELISA was performed using serial two-fold dilutions (from 1:1,600 to 1:102,400) of a serum sample from a human, a horse, and a dog with pythiosis, and a healthy individual (control). The serum dilutions and the obtained EVs of each animal species showed some degree of a linear relationship (i.e., the higher dilution, the lower EV) (**Figure 2**). The result suggests that the protein A/G-based ELISA can quantitate an anti-*P. insidiosum* antibody level for monitoring the clinical course of pythiosis.

All serum samples were also tested in duplicate by the established human-specific ELISA employing goat anti-human IgG antibody conjugated with peroxidase (Bio-Rad Laboratories, California, USA) [15, 23]. The anti-human IgG antibody-based ELISA with an optimal cutoff point (mean EV+4SDs, 15.4; **Table 1; Figure 1B**) did not detect the anti-*P. insidiosum* antibodies in any control sera tested (EV range, 0.0-9.3). This assay can detect the anti-*P. insidiosum* antibodies only in the sera from humans (EV range, 19.8-72.0) and a dog (at a relatively-low EV, 17.9), but not horses (EV range, 6.1-6.6), with pythiosis. The result indicates that the goat anti-human IgG antibody could strongly bind human antibodies and slightly cross-react with dog antibodies. Based on the same serum samples from humans and animals, the anti-human IgG antibody-based ELISA showed equivalent detection specificity (100%), but markedly-lower detection sensitivity (84% vs. 100%), when compared with the protein A/G-based ELISA (**Table 1**).

With a proper cutoff point, all pythiosis sera can ultimately be distinguished from the control sera by the protein A/G-based ELISA (**Table 1; Figure 1A**). Unlike the species-specific ELISAs established here (**Table 1; Figure 1B**) and by other investigators [10–15], the protein A/G-based ELISA is capable of detecting not only the anti-*P. insidiosum* antibodies in a single animal species, but also in humans and various animals (i.e., horses and dogs) those are the most affected hosts of pythiosis. A key feature of our ELISA is the use of protein A/G (rather than a species-specific IgG) as a conjugate, which can bind the anti-*P. insidiosum* immunoglobulins from humans and other animal species. Recently, Intaramat et al. have developed the protein A/G-based ICT, as a rapid and user-friendly assay for diagnosis of pythiosis in humans and animals [16]. They performed a side-by-side performance comparison of their ICT [16] and the species-specific ELISA [11, 13, 15] and showed that ICT has equivalent specificity but relatively-lower sensitivity. One major obstacle for clinical use of ICT is the assay production, which is very complicated and requires a skilled developer and some special reagents that are not available in general clinical laboratories. As opposed to ICT, although ELISA is a multi-step assay and has a longer turnaround time, the production of ELISA is more feasible and requires widely-available reagents. In conclusion, protein A/G-based ELISA has been successfully developed for immunodiagnosis of both human and animal pythiosis. It is a versatile, feasible-to-develop, and functional immunodiagnostic assay.

Limitations

The newly-developed protein A/G-based ELISA overcame the limitation of the previously-established ELISAs [10–15], as it can detect the serum anti-*P. insidiosum* antibodies in multiple host species, including the most affected hosts: humans, horses, and dogs. Other investigators could adopt our protocol to develop their own in-house ELISA. However, there are several limitations in association with the development and performance evaluation of the assay that need to be addressed:

1. The *insidiosum* strain Pi-S was used to prepare the coating antigen. Although this strain has been studied extensively [21, 22, 27], it is not available in the ATCC and CBS-KNAW culture collection centers. Alternatively, any other *P. insidiosum* strains could be used to prepare the coating antigen as well.
2. In the evaluation of our protein A/G-based ELISA, only pythiosis serum samples from 4 horses and a dog were tested, due to the low prevalence of the disease in animals in Thailand.
3. The current study focuses on reporting an optimized protocol for the development of protein A/G-based ELISA. Thus, the diagnostic performance (i.e., detection sensitivity and specificity) of the assay was not extensively investigated, using a variety of serum samples from humans and animals with pythiosis and other diseases.
4. Like the previously-established ELISAs, one extra step to perform the assay is to make a serum dilution (our assay requires 1:1,600 dilution). This step is necessary for reducing the background ELISA signals from the bindings of non-specific antibodies and *insidiosum*'s crude protein extract. In the future, if a purified *P. insidiosum*-specific protein is available for coating the ELISA plate, then an undiluted serum could be directly tested.
5. Pythiosis patients usually present with a chronic infection, as the mean duration of symptoms before receiving medical care is several months [3]. Direct detection of *insidiosum* is an indicator of an active infection. Isolation of *P. insidiosum* from a clinical specimen is difficult and challenging. Antibody detection is an alternative method, providing indirect diagnostic evidence for the *P. insidiosum* infection. Like the other established immunodiagnostic tests for pythiosis, the protein A/G-based ELISA cannot selectively detect an acute-phase infectious marker, such as IgM. Therefore, it cannot directly differentiate between acute, chronic, and previous *P. insidiosum* infections. However, in a patient with the typical features of pythiosis, an elevated or high level of the serum anti-*P. insidiosum* antibodies detected by ELISA suggests an active infection.

Abbreviations

CFA	Culture filtrate antigen
ELISA	Enzyme-linked immunosorbent assay
EV	ELISA value
HA	Hemagglutination

ICT	Immunochromatographic test
ID	Immunodiffusion
OD	Optical density
PBS	Phosphate-buffered saline
SD	Standard deviations

Declarations

Ethics approval and consent to participate

This work was approved by the Committee for Research, Faculty of Medicine Ramathibodi Hospital, Mahidol University (approval number: MURA2019/584).

Availability of data and material

The *P. insidiosum* strain Pi-S (for crude protein preparation) is available upon request.

Funding

The financial supports were obtained from Thailand Research Fund (T. Krajaejun; Grant number: RSA6280092), Faculty of Medicine, Ramathibodi Hospital, Mahidol University (T. Krajaejun; Grant number: CF_61007), and Royal Golden Jubilee Ph.D. Scholarship Program, Thailand Research Fund (C. Jaturapaktrarak; Grant number: PHD/0059/2561).

Acknowledgements

None.

Consent for publication

Not applicable.

Competing interests

None.

Author's contributions

Study design (CJ, TK); Experiment & Methodology (CJ, PP, TL, YK, AL, WP, CY, WT, TK); Data collection (CJ, TL, YK, AL, WP, CY, WT, TK); Data analysis (CJ, PP, TK); Manuscript writing (CJ, TK).

References

1. Gaastra W, Lipman LJ, De Cock AW, Exel TK, Pegge RB, Scheurwater J, et al. *Pythium insidiosum*: an overview. *Vet Microbiol*. 2010;146:1–16.
2. Mendoza L, Ajello L, McGinnis MR. Infection caused by the Oomycetous pathogen *Pythium insidiosum*. *J Mycol Med*. 1996;6:151–64.
3. Krajaejun T, Sathapatayavongs B, Prachartam R, Nitiyanant P, Leelachaikul P, Wanachiwanawin W, et al. Clinical and epidemiological analyses of human pythiosis in Thailand. *Clin Infect Dis*. 2006;43:569–76.
4. Krajaejun T, Lohnoo T, Jittorntam P, Srimongkol A, Kumsang Y, Yingyong W, et al. Assessment of matrix-assisted laser desorption ionization-time of flight mass spectrometry for identification and biotyping of the pathogenic oomycete *Pythium insidiosum*. *Int J Infect Dis IJID Off Publ Int Soc Infect Dis*. 2018;77:61–7.
5. Keeratijarut A, Karnsombut P, Aroonroch R, Srimuang S, Sangruchi T, Sansopha L, et al. Evaluation of an in-house immunoperoxidase staining assay for histodiagnosis of human pythiosis. *Southeast Asian J Trop Med Public Health*. 2009;40:1298–305.
6. Rujirawat T, Sridapan T, Lohnoo T, Yingyong W, Kumsang Y, Sae-Chew P, et al. Single nucleotide polymorphism-based multiplex PCR for identification and genotyping of the oomycete *Pythium insidiosum* from humans, animals and the environment. *Infect Genet Evol*. 2017;54:429–36.
7. Grooters AM, Gee MK. Development of a nested polymerase chain reaction assay for the detection and identification of *Pythium insidiosum*. *J Vet Intern Med*. 2002;16:147–52.
8. Keeratijarut A, Lohnoo T, Yingyong W, Rujirawat T, Srichunrusami C, Onpeaw P, et al. Detection of the oomycete *Pythium insidiosum* by real-time PCR targeting the gene coding for exo-1,3- β -glucanase. *J Med Microbiol*. 2015;64:971–7.
9. Vanittanakom N, Supabandhu J, Khamwan C, Praparattanapan J, Thirach S, Prasertwitayakij N, et al. Identification of emerging human-pathogenic *Pythium insidiosum* by serological and molecular assay-based methods. *J Clin Microbiol*. 2004;42:3970–4.
10. Mendoza L, Kaufman L, Mandy W, Glass R. Serodiagnosis of human and animal pythiosis using an enzyme-linked immunosorbent assay. *Clin Diagn Lab Immunol*. 1997;4:715–8.
11. Grooters AM, Leise BS, Lopez MK, Gee MK, O'Reilly KL. Development and evaluation of an enzyme-linked immunosorbent assay for the serodiagnosis of pythiosis in dogs. *J Vet Intern Med*. 2002;16:142–6.
12. Keeratijarut A, Lohnoo T, Yingyong W, Sriwanichrak K, Krajaejun T. A peptide ELISA to detect antibodies against *Pythium insidiosum* based on predicted antigenic determinants of exo-1,3-beta-glucanase. *Southeast Asian J Trop Med Public Health*. 2013;44:672–80.
13. Santurio JM, Leal AT, Leal ABM, Alves SH, Lübeck I, Griebeler J, et al. Teste de ELISA indireto para o diagnóstico sorológico de pitiose. *Pesqui Veterinária Bras*. 2006;26:47–50.
14. Krajaejun T, Kunakorn M, Niemhom S, Chongtrakool P, Prachartam R. Development and evaluation of an in-house enzyme-linked immunosorbent assay for early diagnosis and monitoring of human pythiosis. *Clin Diagn Lab Immunol*. 2002;9:378–82.

15. Chareonsirisuthigul T, Khositnithikul R, Intaramat A, Inkomlue R, Sriwanichrak K, Piromsontikorn S, et al. Performance comparison of immunodiffusion, enzyme-linked immunosorbent assay, immunochromatography and hemagglutination for serodiagnosis of human pythiosis. *Diagn Microbiol Infect Dis.* 2013;76:42–5.
16. Intaramat A, Sornprachum T, Chanrathonkul B, Chaisuriya P, Lohnoo T, Yingyong W, et al. Protein A/G-based immunochromatographic test for serodiagnosis of pythiosis in human and animal subjects from Asia and Americas. *Med Mycol.* 2016;54:641–7.
17. Jindayok T, Piromsontikorn S, Srimuang S, Khupulsup K, Krajaejun T. Hemagglutination test for rapid serodiagnosis of human pythiosis. *Clin Vaccine Immunol.* 2009;16:1047–51.
18. Prachartam R, Changtrakool P, Sathapatayavongs B, Jayanetra P, Ajello L. Immunodiffusion test for diagnosis and monitoring of human pythiosis insidiosum. *J Clin Microbiol.* 1991;29:2661–2.
19. Mendoza L, Kaufman L, Standard PG. Immunodiffusion test for diagnosing and monitoring pythiosis in horses. *J Clin Microbiol.* 1986;23:813–6.
20. Krajaejun T, Imkhieo S, Intaramat A, Ratanabanangkoon K. Development of an immunochromatographic test for rapid serodiagnosis of human pythiosis. *Clin Vaccine Immunol.* 2009;16:506–9.
21. Rujirawat T, Patumcharoenpol P, Lohnoo T, Yingyong W, Lerksuthirat T, Tangphatsornruang S, et al. Draft Genome Sequence of the Pathogenic Oomycete *Pythium insidiosum* Strain Pi-S, Isolated from a Patient with Pythiosis. *Genome Announc.* 2015;3:e00574-15.
22. Rujirawat T, Patumcharoenpol P, Lohnoo T, Yingyong W, Kumsang Y, Payattikul P, et al. Probing the Phylogenomics and Putative Pathogenicity Genes of *Pythium insidiosum* by Oomycete Genome Analyses. *Sci Rep.* 2018;8:4135.
23. Lohnoo T, Yingyong W, Kumsang Y, Payattikul P, Jaturapaktrarak C, Chailurkit L-O, et al. Seroprevalence of anti-*Pythium insidiosum* antibodies in the Thai population. *Med Mycol.* 2019;57:284–90.
24. Chaiprasert A, Samerpitak K, Wanachiwanawin W, Thasnakorn P. Induction of zoospore formation in Thai isolates of *Pythium insidiosum*. *Mycoses.* 1990;33:317–23.
25. Lohnoo T, Jongruja N, Rujirawat T, Yingyong W, Lerksuthirat T, Nampoon U, et al. Efficiency comparison of three methods for extracting genomic DNA of the pathogenic oomycete *Pythium insidiosum*. *J Med Assoc Thai.* 2014;97:342–8.
26. Badenoch PR, Coster DJ, Wetherall BL, Brettig HT, Rozenblds MA, Drenth A, et al. *Pythium insidiosum* keratitis confirmed by DNA sequence analysis. *Br J Ophthalmol.* 2001;85:502–3.
27. Krajaejun T, Rujirawat T, Kanpanleuk T, Santanirand P, Lohnoo T, Yingyong W, et al. Biochemical and genetic analyses of the oomycete *Pythium insidiosum* provide new insights into clinical identification and urease-based evolution of metabolism-related traits. *PeerJ.* 2018;6:e4821.

Table

Table 1. Correlation between the cutoff points and the diagnostic performance (i.e., sensitivity, specificity, and accuracy) of the protein A/G-based and anti-human IgG antibody-based ELISAs for detection of the anti-*P. insidiosum* antibodies in serum samples.

Cutoff point ^a	Protein A/G-based ELISA				Anti-human IgG antibody-based ELISA			
	Cutoff value (EV) ^b	Sensitivity ^c	Specificity ^d	Accuracy ^e	Cutoff value (EV)	Sensitivity	Specificity	Accuracy
Mean EV + 5SDs	13.4 ^f	100.0	100.0	100.0	18.6	80.0	100.0	93.3
Mean EV + 4SDs	11.3	100.0	98.0	98.7	15.4 ^g	84.0	100.0	94.7
Mean EV + 3SDs	9.3	100.0	98.0	98.7	12.2	84.0	100.0	94.7
Mean EV + 2SDs	7.2	100.0	96.0	97.3	9.0	84.0	98.0	93.3

Footnotes

^a Cutoff point = mean EV of the control sera plus 2, 3, 4 or 5 SDs

^b ELISA value (EV) = OD value of a serum sample / OD value of the negative control serum

^c Sensitivity (%) = true positive samples / (true positive + false negative samples) x 100

^d Specificity (%) = true negative samples / (true negative + false positive samples) x 100

^e Accuracy (%) = (true positive + true negative samples) / all samples x 100

^f Optimal cutoff point of the protein A/G-based ELISA

^g Optimal cutoff point of the anti-human IgG antibody-based ELISA

Figures





Figure 1

Dot plot of the ELISA values (EV) of pythiosis (n = 25) and control (n = 50) serum samples examined by (A) the protein A/G-based and (B) the anti-human IgG antibody-based ELISAs. Different colors indicate the EVs of sera from different animal species. Dashed lines represent the optimal ELISA cutoff points (see Table 1).




Figure 2

ELISA values (Y-axis) of serial two-fold dilutions (from 1:1,600 to 1:102,400; X-axis) of a serum sample from a human (blue), a horse (orange), and a dog (gray) with pythiosis, and a healthy individual (served as a control; yellow). The working serum dilution for protein A/G-based ELISA testing is 1:1,600 (see text for details).